



PennWhite

Foamdoctor F2875

foam control agent for fermentation processes

technical data sheet

foam control agents - food grade

Foamdoctor F2875 is a concentrated foam control agent for use in various fermentation applications where foaming is an inherent problem. **Foamdoctor F2875** has shown to be effective in the control of foam produced in wide range of fermentation applications, including citric acid, yeast, antibiotic, alcohol, glutamic acid, enzyme, protein, pharmaceutical and nucleotide production.

TYPICAL PROPERTIES

Appearance:	Colourless liquid.
Specific Gravity:	Ca 1.0
Active Content:	100%.
Viscosity:	Ca 1000 cps @ 25°C.
Dispersibility:	Dispersible in cold water (15°C).

APPLICATION

Foamdoctor F2875 is of low order toxicity, has a neutral odour and taste and has relevant food grade approvals.

Foamdoctor F2875 can be safely utilised in fermentation processes where ultra filters are used since it is dispersible/self-emulsifiable in water. It is most soluble on cooling to below 20C, hence virtually all traces of de-foamer can be removed by washing with cold water.

Foamdoctor F2875 can be used in extremely small amounts and combined with the potential for increased productivity makes **Foamdoctor F2875** an extremely cost effective processing aid.

Foamdoctor F2875 is silicone and silica free and therefore does not create filtration or oxygen transfer problems that can lead to reduced growth rates. Such problems can be created with products containing silicones.

Foamdoctor F2875 can be easily sterilised by normal heating and/or autoclaving methods without affecting product stability or performance.

Foamdoctor F2875 has been scientifically designed to optimise the use of polyalkylene glycol technology, in order to meet the requirements of the wide range of applications in the fermentation industry

For most effective use of **Foamdoctor F2875** it should be added continuously via a metered pump. The dosage level of **Foamdoctor F2875** required will normally be in the region of 10 –50 ppm of the volume of culture medium. The exact dosage required will depend on the foaming propensity of the culture medium.

Foamdoctor F2875 is designed to comply with current UK (MAFF) regulations as well as one or more of FDA CFR 21 Nos 173.340, 175.105, 175.300, 176.170, 176.180, 176.200, 177.1220, 177.2260; BGA Chapters 14 and 36 and the French "Repression des Fraudes".

PACKAGING

Foamdoctor F2875 is available in 25 litre kegs, 205 litre drums, 1000 litre IBC containers and bulk tanker.

STORAGE & SHELF LIFE

Foamdoctor F2875 should ideally be stored for use in temperature conditions between minus10 °C and 40°C. **Foamdoctor F2875** has a shelf life of at least twelve months when stored in the recommended conditions in original unopened containers.

(OPA2/F2)

TDS: 101

Issue: 04

Revised: 06/03/2014

Page 1 of 1

For additional technical or safety advice call us
on 01606 734820 or email info@pennwhite.co.uk

PennWhite Ltd
Aston Way
Midpoint 18 Business Park
Middlewich, CW10 0HS



Telephone: +44 (0)1606 734820
Fax: +44 (0)1606 837867
Email: info@pennwhite.co.uk
Website: www.pennwhite.co.uk